
***Electron Energy-Loss Spectroscopy (EELS)
for the Hitachi HD-2000***


Outline

- Introduction to Electron Energy-Loss Spectroscopy (EELS)
- EELS Spectrum Imaging
- EELS on the HD-2000

Introduction to EELS

- ❑ How EELS works
- ❑ What can EELS do?
- ❑ Spectrum information
- ❑ Quantification
- ❑ Comparisons to X-ray microanalysis

How EELS works



- (S)TEM probe electrons travel through a thin specimen
- Probe electrons lose energy due to their interaction with the specimen (*inelastic scattering*)
- The energy-losses are characteristic of the elements and chemistry of the specimen
- An EELS Spectrometer can disperse the probe electron beam according to its lost energies into a spectrum

What can EELS do?


- EELS can provide information about:
 - Elemental composition of all elements $Z > 1$
 - Chemical composition
 - Specimen thickness

- The high sensitivity of EELS makes it ideal for mapping


- Further capabilities:
 - *Chemical bonding*
 - *Band structure*
 - *Valence and conduction electron density*
 - *Atom-specific radial distribution of near neighbors*
 - *Polarization/dielectric response*
 - *$Z=1$ quantification*

Basic Physics of EELS

- (S)TEM probe electron donates energy to core or valence electrons in the sample
- EELS analyzes the energy lost in the probe electrons



EELS Spectrum Features




EELS Spectrum Quantification

- ❑ Calculate edge intensities by removing background below edge
- ❑ Convert edge intensities into concentrations using cross-sections (σ)
- ❑ No standards needed!

- ❑ Absolute quant
 - Use edge and low-loss intensity
 - Result is in atoms or atomic density

- ❑ Relative quant
 - Uses ratio of edge intensities
 - Result is in atomic fractions



Automated EELS elemental analysis


Gatan software performs fully automated absolute and relative elemental quantification

Spectrum Display

Atomic ratios



Ba/O: 0.16 ± 0.024 Cu/O: 0.10 ± 0.017

Cursor: 1048.2 Width: 0.48422 Counts: 0



Chemical analysis

- Changes in chemical bonding produce changes in the shape of EELS edges
 - Fingerprinting is the technique of identifying a chemical structure from the edge shape
 - Quantification is possible with references




EELS Resolution

- EELS information has excellent energy and spatial resolution
 - Spatial resolution for EELS edges is limited only by probe size
 - Energy resolution is limited by the probe energy distribution

- Typical energy resolution (*depends on gun current*)
 - LaB₆ : >1 eV
 - Warm FEG: 0.6-0.9 eV
 - Cold FEG: 0.35-0.5 eV

- Cold FEG STEM has best energy and spatial resolution



EDXS and EELS comparisons

- EELS and X-ray are complementary techniques
 - EELS excels at lower Z
 - X-ray excels at higher Z
 - EELS and X-ray data can be acquired simultaneously

- Quantification
 - EELS processing doesn't need standards (no k-factors)
 - Standards-based processing is easier with X-ray (currently less developed for EELS since it is less widely used)
 - EELS has greater sensitivity for $Z < 30$ (and most of periodic table)
 - EELS can do absolute quant or relative quant

- Mapping
 - EELS mapping is much faster because of greater signal generation and detection efficiency
 - EELS maps have greater S/N and more pixels for the same mapping times

EDXS and EELS comparisons

(S)TEM EDS X-ray

- ❑ X-rays provide elemental information only
- ❑ Inefficient signal collection; inefficient low Z signal generation & detection
 - Slow mapping or poor S/N
- ❑ X-ray spectra can contain artifact information from column and other parts of sample
- ❑ High detection efficiency for higher Z elements
 - Poor sensitivity to $Z < 10$
- ❑ Energy resolution $> 120\text{eV}$ causes frequent overlaps
- ❑ No sample thickness limitations

(S)TEM EELS

- ❑ EELS provides elemental, chemical & dielectric information
- ❑ Very efficient in all aspects
 - Higher sensitivity to most elements
 - Very fast mapping technique
- ❑ EELS information is highly localized and does not contain sample or column artifacts
- ❑ High detection efficiency for lower Z elements
 - Poor sensitivity to a few high Z elements
- ❑ Energy resolution $0.3\text{-}2\text{eV}$ gives far fewer overlaps (overlaps when edges $\sim < 30\text{eV}$ apart)
- ❑ Sample thickness is important - should be less than $\sim 100\text{nm}$ @ 200keV

Spectrum imaging

- ❑ What is Spectrum Imaging?

- ❑ Gatan spectrum imaging
 - Overview
 - Correcting for instabilities
 - Basic processing

- ❑ Examples

What is Spectrum Imaging?


- Spectrum imaging is the technique of acquisition, storage, processing, and analysis of spectroscopic data at each pixel in a digital image
- Conventional image pixels contain one intensity value only
- Spectrum image pixels contain complete spectra
- Produces quantitative elemental maps and profiles

Spectrum image schematic:

$\Delta x, \Delta y$ Spatial dimensions
 ΔE Energy-loss dimension


Vertical columns are spectra

Horizontal planes are energy filtered images



Acquisition of Spectrum Images

- A spectrum image is typically acquired in STEM mode by stepping a focused electron probe from one pixel to the next
- The spectrum image data cube is filled one spectrum column at a time
- In STEM it is possible to collect EELS, X-ray or both spectra simultaneously
- Use of the DF or SE signal during acquisition permits spatial drift correction and assurance that information is coming from the desired area




Gatan Spectrum Imaging Overview

- Powerful and highly automated software for producing quantitative maps and line profiles
 - Supplied as plug-ins to Gatan DigitalMicrograph scientific-grade image processing software
 - Intuitive acquisition and analysis features
 - EELS spectrum imaging can map up to 50 pixels/second
 - Corrects for detector artifacts and system instabilities
 - X-ray spectrum imaging available in Fall '99 for many brands

- Requires Gatan DigiScan hardware
 - Controls beam raster
 - ◆ 16 bit DACs (max image size 8K*8K)
 - Simultaneous acquisition of multiple signals (BF, DF, SE, etc.)
 - ◆ 12-bit ADC samples at 5MHz, can oversample to >16 bits
 - Flexible control of all acquisition parameters (size, speed, aspect ratio ...)

Correcting system instabilities

- Gatan software can automatically correct:
 - All detector artifacts
 - System instabilities
 - ◆ Spatial drift
 - ◆ Energy drift
 - ◆ Current drift
 - User needs only to enable correction
 - User can modify defaults for improved control
 - ◆ Reduce dose to analyzed area
 - ◆ Increase correction frequency



Basic spectrum image processing










Image Status
Image B
Type: Real 4
Size: 64 x 90 x 1019
X: 46.8 nm
Y: 65 nm
Zc: 281 eV
Value: 2.67e+06 e-
Slice
Low: 276 eV
High: 287 eV
Width: 10.8 eV
Display Center





Example - Spectrum imaging




- ❑ EELS spectrum image [160,90] pixels, 1022 energy channels
- ❑ Acquired on HD-2000 at 200keV

Example - Line profiles

- 100 points acquired in < 1 min
- Acquired on HD-2000 at 200keV



Example - Spectrum imaging



EELS on the HD-2000

- ❑ Main DigiPEELS features
- ❑ DigiPEELS on the HD-2000
- ❑ Configurations for HD-2000
- ❑ Major selling points of EELS on HD-2000

Main DigiPEELS features

- ❑ DigiPEELS is a parallel-detection energy-loss spectrometer
 - EELS spectrum channels acquired in parallel


- ❑ All spectrometer features fully computer controlled

- ❑ Linear photodiode array detector
 - 1024 channels
 - 16-bit digitization
 - Maximum of 50 spectra/sec

- ❑ Automated EELS processing software



DigiPEELS on HD-2000



Not to scale

EELS is Simplified on HD-2000

- ❑ The most difficult issues of EELS data acquisition have been simplified or eliminated on HD-2000 systems !!!
 - Due primarily to lack of TEM viewing chamber and short camera lengths
 - Virtually all scattering angles acquired which simplifies quantification

Tips for good EELS data acquisition

- ❑ Work in diffraction or STEM modes
 - Stay away from image mode - chromatic aberration risks
- ❑ Control your collection angles
 - Know your collection and convergence angles
 - Don't use an aperture too small to admit desired electrons
 - Center diffraction on spectrometer aperture
- ❑ Use valid detector artifact corrections

Not necessary on HD-2000 + DigiPEELS systems!

DigiPEELS use on HD-2000 Details

- ❑ On HD-2000, EELS is done with projector lens turned off
 - Camera length is extremely short
 - Spectrometer focussing virtually unnecessary
 - Switching between EELS apertures unnecessary
 - Essentially all of the inelastic signal is collected ($>100\text{mrad}$)
 - Precise knowledge of collected scattering angles unnecessary
 - Precise centering of beam in DigiPEELS unnecessary

- ❑ HD-2000 has no projector chamber
 - Very good immunity to external fields that affect EELS resolution

- ❑ Microscope alignment
 - Installation engineer mechanically aligns DF detector to DigiPEELS
 - Thereafter, user alignment of microscope for imaging also aligns DigiPEELS

EELS Configurations for HD-2000

- DigiPEELS (*available on HD-2000 in Fall '99*)
 - HD-2000 configured DigiPEELS
 - HD-2000 dark field detector with mechanical alignments and hole $\geq 3\text{mm}$
 - Recent Macintosh or Windows (NT/98) computer
 - Options
 - ◆ Diffraction pattern observation camera mounted above PEELS
 - TV rate camera - Gatan model 692
 - 1K*1K large pixel cameras: Gatan models 792/794

- EELS Spectrum Imaging (*available on HD-2000 in Fall '99*)
 - HD-2000 configured DigiPEELS
 - HD-2000 configured DigiScan (controls beam scan & ultra high-quality imaging)
 - STEM EELS spectrum imaging & DigitalMicrograph image processing software
 - Options
 - ◆ X-ray spectrum imaging software (not all brands available Fall '99)

Major selling points of EELS on HD-2000

- Point Analysis
 - Low-Z elements that X-ray does poorly (Li,Be,B,C,N,O,F)
 - Detection of very small concentrations (<100 ppm possible for most elements)
 - Use when X-ray sample/column/detector artifacts are a problem
 - ◆ EELS can be used on FIB made specimens without the artifacts X-ray has from around FIB box
 - Use when overlaps are a serious problem (such as separating N, Ti, O)

- Mapping with spectrum imaging
 - All advantages of point analysis
 - Very fast mapping times

- Simplicity
 - HD-2000 has made EELS more simple and accurate than ever

Contacts at Gatan

- Marketing / Sales Assistance for HD-2000 DigiPEELS & Spectrum Imaging
 - Dr. Gerald Kothleitner
EELS Product Manager
Gatan R&D
5933 Coronado Lane
Pleasanton, CA 94588

(925) 463-0200
(925) 463-0204 {FAX}
geraldk@gatan.com